

Isolator pressure alarm

SAFETY ALARM SYSTEM



Process Containment

Keeps operators and isolators safe

Our isolator pressure alarm reads and displays the pressure inside your isolator, alerting you via both a sound and light alarm if the pressure is either too high or has dropped too low.

Designed to integrate with any type of flexible isolator, the pressure alarm can be used alongside ONFAB's range of fan control systems, helping to secure the integrity of the isolator and protect users.

The alarm is simple to use, with a clear read out and is also lightweight, making it easy to place and move.



Isolator pressure alarm monitors pressure levels in your isolator

KEY FEATURES

- Controls both under and over pressure
- Adjustable to your processing requirements
- Clear pressure read out
- Alarm uses both sound and light to draw attention
- Lightweight
- Easy to use
- Integrates with any flexible isolator



Raises alarm if pressure is too high or too low

QUALITY STATEMENT

ONFAB was founded in 2004 and has grown to become a global leader in the bespoke design, manufacture and installation of flexible containment systems for the pharmaceutical and healthcare industries. Our team of engineers are committed to the highest standards of excellence, meeting the requirements of ISO 9001 ISOQAR. From our production facilities in the UK and Europe, we advise on protection for your workforce and how to streamline and validate the manufacturing process.



FS 753783

4905/02/0322



Onfab

6-7 Prosperity Court | Middlewich, Cheshire CW10 0GD
+44 (0) 1606 832080 | enquiries@onfab.co.uk | www.onfab.co.uk