

**ONFAB****Process Containment**

# Pop-up isolator

**LAB CONTAINMENT**

## Fast, safe handling for high potent pharmaceuticals

This instant pop-up benchtop isolator transforms your lab into a safe, enclosed working environment for handling high potent pharmaceuticals. With quick and easy set up, it's ideal for small scale powder, liquid, instrument and inert-gas operations. It can be used for a variety of applications including:

- IPC (In Process Control)
- Weighing and dispensing in a controlled environment
- Material QA testing
- Instrument protection and housing
- Granulation
- Solution preparation



*Pop-up safety isolator*

### KEY FEATURES

- Set up in less than 2 minutes
- Typical size 770mm long x 620mm wide x 740mm high
- 2 x glove port access
- Large gas-prof zip for material/equipment introduction and removal
- Wide choice of anti-static materials including FDA approved films
- Integrated P3 filter for environmental and operator protection
- Integrated service sleeve for utility cables
- Easy to clean with IPA wipes
- Ergonomic design
- Option to buy separate stainless steel or PTFE isolator base

4487/02\_0322

**ONFAB****Onfab**

6-7 Prosperity Court | Middlewich, Cheshire CW10 0GD  
+44 (0) 1606 832080 | [enquiries@onfab.co.uk](mailto:enquiries@onfab.co.uk) | [www.onfab.co.uk](http://www.onfab.co.uk)

# Pop-up isolator

LAB CONTAINMENT



**ONFAB**


Process Containment



*The pop-up isolator can be set up in less than 2 minutes*



*Features FDA approved films*



*Safe handling of high potent pharmaceuticals*



*Suitable for various lab based containment applications*



FS 753783

## QUALITY STATEMENT

ONFAB was founded in 2004 and has grown to become a global leader in the bespoke design, manufacture and installation of flexible containment systems for the pharmaceutical and healthcare industries. Our team of engineers are committed to the highest standards of excellence, meeting the requirements of ISO 9001 ISOQAR. From our production facilities in the UK and Europe, we advise on protection for your workforce and how to streamline and validate the manufacturing process.

4487/02\_0322



**ONFAB**

Onfab

6-7 Prosperity Court | Middlewich, Cheshire CW10 0GD

+44 (0) 1606 832080 | [enquiries@onfab.co.uk](mailto:enquiries@onfab.co.uk) | [www.onfab.co.uk](http://www.onfab.co.uk)