

CLEANROOM PPE

SUITE® GOGGLES



DESCRIPTION

The Suite® Sterile Goggles are a revolutionary product for the sterile environment. These goggles are single-use to prevent cross contamination at a very cost-effective price. Individually bagged in clean packaging and gamma irradiated to ensure sterility. The gamma irradiated Goggles have an exceptional anti-fog lens along with micro-breathable top vents utilizing direct or indirect ventilation for breathable comfort without fogging. They offer see-through peripheral vision, provide a secure comfortable fit and pass all ANSI-Z87.1 standards for safety.

PRODUCT INFORMATION

Part Numbers:

1. A-50G0G-IV-S
2. A-50B-IR

Packaging:

1. 1 per heat-sealed pack, 50 packs per double case lined carton
2. 1 per heat-sealed pack, 50 packs per double case lined carton
50 packs per double case lined carton

Description:

1. Anti-fog, validated sterile, single-use goggle with indirect vents
2. Anti-fog validated sterile, single-use goggle with direct ventilation

Sterility: Gamma Irradiated



PRODUCT PERFORMANCE

Type	Method	Result
Safety	ANSI-Z87.1	Passed
Sterility	AAMI TIR33:2005	Validated Sterile to a 10 ⁻⁶ Sterility Assurance Level

BP-STE-0008-A

