

CLEANROOM WIPERS

SUITE® SC12 WIPERS



DESCRIPTION

Suite® wipers are engineered for high-performance. Uniquely manufactured from synthetic high-density Ultra Fibers, Suite® wipers are non-particulating, snag resistant, and offer unrivaled absorbency. They deliver maximum performance and durability when used as part of your cleaning and/or disinfecting protocols.

- Double welded seam to reduce particle generation in comparison to traditional knife cut or single welded seam.



PRODUCT INFORMATION

Part Numbers:

- **Non-Sterile:** W-SC12-25
- **Sterile:** W-SC12-5S

Size: 12" x 12"

Color: White

Packaging:

- **Non-Sterile:** 200 wipers/case (25 wipers/bag, 8 bags/case)
- **Sterile:** 200 wipers/case (5 wipers/bag, 4 bags/pack, 10 packs/case)



PRODUCT PERFORMANCE

| Type | Method | Result |
|-----------|-----------------|---|
| Particles | IEST-RP-CC03 | Acceptable for ISO Class 5 Use |
| Sterility | AAMI TIR33:2005 | Validated Sterile to a 10 ⁻⁶ Sterility Assurance Level |

Reduction and Transfer: Wiping with Suite® material has been demonstrated to *effectively remove surface contamination. In addition, the carry-over (cross contamination potential) after wiping with the Suite® material is minimized as contaminants are fully retained by the material.

*Lab testing available upon request.

BP-STE-0008-A

