


Custom Sterile Zip Lock Bags

- Bags are validated sterile to a 10^{-6} S.AL. in accordance with AAMI guidelines. Quarterly bioburden and dose audits are performed to assure sterility.
- Manufactured from high clarity low-density polyethylene in an ISO Class 5 cleanroom.
- Size, gusset and mil thickness can be made to your specifications.
- Available in stock gauge of 2.5 mil in sizes 4"x6", 6"x8", 8"x10", 10"x12", 12"x15", and 18"x24"
- Bags available in a Zip Clean NY COEX clear material.
- Stocked sizes are available for immediate shipment.
- All final product is double bag packaged to further assure product cleanliness and sterility.
- All final goods are identified by a lot number that allows for traceability back through manufacturing and sterilization.
- Gamma irradiated industrial zip lock bags also available at a lower cost.



Product Performance

Test	Method	Result
Sterility	ANSI/AAMI/ISO 11137	Validated Sterile to a 10^{-6} Sterility Assurance Level