

CLEANROOM SUPPLIES

SUITE® STERILE TRASH BAGS AND STAND



DESCRIPTION

Validated sterile polyethylene trash bags are available in three stocked sizes. They are manufactured from high-clarity polyethylene in an ISO Class 5 cleanroom and are validated sterile to a sterility assurance level of 10^{-6} . Each bag is double-bagged to ensure aseptic entry into sterile core areas. Custom sizes are also available to meet specific sterile bag requirements. These bags pair perfectly with our stainless steel X-stand, which is autoclavable and made from chemical-resistant 304 stainless steel. Both the bags and X-stand are proudly made in the USA.

PRODUCT INFORMATION

Part Numbers:

- **Autoclavable X-Stand:** E-SS3340ST
- **Sterile Trash Bags for X-Stand:**
 1. PB3340-4-IR: 33x40, 4 ml
 2. B-3648-1S: 36x48, 1ml
- **Additional Sterile Polyethylene Trash Bags:**
PB2436-4-S: 24x36, 4ml

Packaging:

- **Autoclavable X-Stand:** 1 stand per box
- **Sterile Trash Bags for X-Stand:**
 1. 25/bag, 4 bags/case
 2. 25/pack, 10 packs/case
- **Additional Sterile Polyethylene Trash Bags:**
25 bags per pack, 4 packs/case



PRODUCT PERFORMANCE

Test	Method	Result
Sterility	ANSI/AAMI/ISO 11137	Validated Sterile to a 10^{-6} Sterility Assurance Level

BP-STE-0015-A

