

# CLEANROOM MOPS

## SUITE® TOUCHLESS REMOTE DUTY PATENT PENDING



### DESCRIPTION

Suite® Touchless Remote Duty mops eliminate operator glove handling of both clean and sterile or dirty pads while also reducing time spent changing pads. The compact design allows users to operate in more confined and restricted spaces such as RABS units, small vent hoods, BSC's, Isolators, and other difficult to reach places. The engineered Suite® fiber provides consistent and liberal application of your solutions in a non-particulating and snag resistant dynamic material.

### PRODUCT INFORMATION

#### Part Numbers:

- **Sterile:**
  1. Y-TLM-511-1S
  2. Y-TLM-511-25S
- **Non-Sterile:** Y-TLM-511-25
- **Handle:** Y-1384-M
- **Frame:** Y-10510-AANK
- **Strips:**
  1. Y-10510-50IR
  2. Y-10510-50

#### Description:

- **Handle:** Ergonomic autoclavable anodized aluminum and plastic with handle that extends from 19" to 35".
- **Frame:** Anodized aluminum autoclavable 10" Suite touch-less frame.

#### Packaging:

- **Sterile:**
  1. 1 mop per bag, 5 bags per pack, 5 packs per case
  2. 5 mops per pack, 25 mops per case
- **Non-Sterile:** 5 mops per pack, 25 mops per case



### SPECIFICATIONS

Type	Method	Result
Particles	IEST-RP-CC03.4   Helmke Drum	Acceptable for ISO Class 5 Use   Category 1
Sterility	AAMI TIR33:2005	Validated Sterile to a 10 <sup>-6</sup> Sterility Assurance Level

**Reduction and Transfer:** Wiping with Suite® material has been demonstrated to \*effectively remove surface contamination. In addition, the carry-over (cross contamination potential) after wiping with the Suite® material is minimized as contaminants are fully retained by the material.

\*Lab testing available upon request.

BP-STE-0008-A

