

PRE-SATURATED POLYPROPYLENE WIPES

SPECIFICATION SHEET



DESCRIPTION

BenchmarkGMP Pre-Saturated Polypropylene Wipes are made from 100% melt blown polypropylene and they are saturated with a 70% isopropyl alcohol and 30% WFI DI water mixture giving them excellent durability and saturation. The pouches are individually bagged within a case liner to ensure sterility while remaining an economical option for both cleaning and disinfecting.



APPLICATIONS

- The material makeup of these wipes means they are low particle and fiber generating while maintaining durability.
- Ideal for use in cleanrooms, on equipment and in other critical manufacturing areas.



SPECIFICATIONS

Type

Particles
Fibers

Method

Orbital Shake Test IEST-4.3-6.1.4
Orbital Shake Test IEST-4.3-6.1.4

Results

$48 \times 10^6/\text{m}^2$
 $7.6 \times 10^4/\text{m}^2$

Part Number	W-BP7030-911
Size	9" x 11"
Packaging	50 wipes/pack, 24 packs/case
Material	100% Meltblown Polypropylene
Solution	70% IPA 30% WFI DI Water
Basis Weight	35 g/m ²
Sorptive Capacity	260 ml/m ²
Sorptive Rate	< 1 second