

# BOXED BLUE NITRILE GLOVES

## SPECIFICATION SHEET



### DESCRIPTION

These durable and high-quality boxed blue nitrile gloves are strong yet comfortable. They feature a micro-textured fingertip and beaded cuff for maximum comfort and control. The gloves are made with material to reduce the risk of allergic reactions.

These boxed gloves are available in sizes Small to X-Large in boxes of 100 gloves – 10 boxes per carton.

### PERFORMANCE

- These gloves meet the requirements of ASTM D6319-10
- Gloves are inspected and tested in accordance with ISO2859-1:1999
- Samples are inspected and tested per ASTM D6319-10



Part Numbers	G-BP10B-S   G-BP10B-M   G-BP10B-L   G-BP10B-XL
Description	Powder Free Nitrile Gloves with Micro-Textured Finish
Color	Blue
Texture	Textured Fingertips
Weight	4.5g $\pm$ 0.2g
Cuff	Beaded
Type	Ambidextrous
Packaging	100 gloves per box, 10 boxes per carton
Length	230mm minimum (all sizes) – (ASTM D6319-10)
Thickness	3.5 mil
AQL Freedom from Holes	1.5 (ASTM D6319-10)
Tensile Strength (Mpa)	14 Before Aging   14 After Aging (ASTM D6319-10)
Elongation %	500 Before Aging   400 After Aging (ASTM D6319-10)

\*Sourced from a Small Minority-Owned Supplier (Hub Zone and SBA8(a))

V.9.BB

