

## PRECISION CLEAN® I

Gauges from .0015" to .010"  
Custom flat bags sizes from 2" x 3" to 47" x 120"  
Custom gusseted sizes (please contact us for capabilities)  
Available printed up to eight colors  
Additional available film and bag variations  
Custom tubing, centerfold and sheeting are available



Precision Clean I® bags are manufactured in our ISO Class 4 Cleanroom, available printed and unprinted, and suitable for gamma and ETO sterilization. The PCI resin and film base have been carefully chosen with cost in mind. PCI is specifically designed to service less critical cleanroom fill operations that need films with Surface Cleanliness Level 200 or higher (Cleanrooms that are less controlled than ISO Class 4). PCI has excellent clarity.