

Return To Controlled Environment Class 100 ISO 5 Gloves Products

CLASS 100 CONTROLLED ENVIRONMENT STERILE LATEX GLOVES

SGPF SERIES

CTI's SGPF Series are Sterile Gloves. Like all our gloves, they are manufactured for outstanding strength, stretch and fit. They offer a textured surface and are double chlorinated reducing contact sensitizers and protein content.

SGPF Series gloves are lot traceable, 12" long and hand-specific. They have been tested to meet and/or exceed all standards and requirements. Available in sizes 6 to 10.



PRODUCT FEATURES

- Textured Surface
- Low Protein and Contact Sensitizers
- Double Chlorinated
- Tapered, Beaded Cuff

- Lot Traceable

- Folded Cuff for Aseptic Donning

- 12" Length, Hand Specific

- Total Protein Content <50ug/g of Glove (ppm) Using Modified Lowry Method

- Leak Tested

- Meets ASTM D3578-05 standards

- Sampled in accordance with MIL STD-105E

- Particle & Extractable Testing in Accordance with IES-RP-CC-005.2

- Qualified for Use in the Pharmaceutical Industry

- Gamma Irradiated Sterilization Minimum Dose 2.5 Mrads

- Sterilized in Accordance with AAMI Standards

SPECIFICATIONS & PACKAGING

Length: 12" (+/-0.5)

Thickness: 6.5 mil (+/-1)

Tensile Strength: 400 psi/23 mpa

Elongation: 745%

Packaging: HDPE/Poly Wallet, 1 Pair/Pouch, 50 Pairs/Box, 200 Pairs/Case

Cat No./Sizes:

SGPF600 – 6.0

SGPF650 – 6.5

SGPF700 – 7.0

SGPF750 – 7.5

SGPF800 – 8.0

SGPF850 – 8.5

SGPF900 – 9.0

SGPF1000 – 10.0