

2.1 125ml PET (100125-05)

2.1.1 Product Specification Sheet 125ml PET

PharmaTainer specifications 125mL PET

125 ml PET PharmaTainer (PharmaPack)			Catalogue number 100125-05		
Bottle part number	100125-single	USP Class VI	ADCF	DMF No.	
Bottle material	Cellon resin code 001	Yes	Yes	15205	
Nominal volume	125ml				
Brim capacity	196ml	Tolerance	Shape (footprint) Square		
Bottle dimensions w/o cap	H95 x W58.8 x D58.8mm	+/- 1.2mm			
Bottle weight w/o cap	50.0g	+/- 3%			
Bottle thread O.D.	44mm	+/- 0.3mm			
Bottle neck I.D.	37mm	+/- 0.3mm	Accuracy		
Graduations	25ml – 125ml	25ml intervals	+/- 5%		Endotoxin
Sterility	C60 γ irradiated X-Ray	SAL 10 ⁻⁶			No detectable Endotoxin Certified by Lot
Bottle decoration	Graduation scale	Lot number	Expiration date		Data matrix
Bottle colour	Clear				
Traceability	Lot number	Data matrix	Serial number		
Closure part number	300048C-99		USP Class VI / EP	ADCF	DMF No.
Closure material	Cellon resin code 003		Yes	Yes	1646
			Tolerance		
Closure weight	12.6g		+/-2g		
Closure dimensions	H 25mm O.D 49mm		+/- 0.3mm		
Closure colour	Natural				
Packaging materials	Carton	Level 3 Carton liner	Level 2	Level 1 Sterile barrier	
	4mm single corrugate, top and base protection sheet	50µ LDPE monolayer film tape sealed bag	100µ LDPE monolayer film Heat sealed bag	140µ PE/PA multilayer film vaccum bag heat sealed	
Carton dimensions	Carton weight		Units per carton	Units per sterile barrier	
L 620 – W340 – H 270 mm	7.86kg		100	5 (PharmaPack 5)	
	Dimensions L x W x H	Weight (inc. pallet)	Cartons per pallet	Units per pallet	
Euro pallet	1200x800x2190mm	207kg	34	2400	
Industry pallet	1200x1000x1820mm	257kg	30	3000	