

2.3 1L PET (101000-05)

2.3.1 Product Specification Sheet 1L PET

PharmaTainer specifications 1L PET

1L PET PharmaTainer (PharmaPack)			Catalogue number 101000-05	
Bottle part number	101000-single	USP Class VI	ADCF	DMF No.
Bottle material	Cellon resin code 001	Yes	Yes	15205
Nominal volume	1000ml			
Brim capacity	1213ml	Tolerance		Shape (footprint)
Bottle dimensions w/o cap	H187 x W101.2 x D101.2 mm	+/- 1.2mm		Square
Bottle weight w/o cap	165.8g	+/- 3%		
Bottle thread O.D.	44mm	+/- 0.3mm		
Bottle neck I.D.	37mm	+/- 0.3mm	Accuracy	
Graduations	100ml – 1000ml	100ml intervals	+/- 5%	Endotoxin
Sterility	C60 γ irradiated X-Ray	SAL 10 ⁻⁶		No detectable Endotoxin Certified by Lot
Bottle decoration	Graduation scale	Lot number	Expiration date	Data matrix
Bottle colour	Clear			
Traceability	Lot number	Data matrix	Serial number	
Closure part number	300048C-99	USP Class VI / EP	ADCF	DMF No.
Closure material	Cellon resin code 003	Yes	Yes	1646
		Tolerance		
Closure weight	12.6g	+/-2g		
Closure dimensions	H 25mm O.D 49mm	+/- 0.3mm		
Closure colour	Natural			
Packaging materials	Carton	Level 3 Carton liner	Level 2	Level 1 Sterile barrier
	4mm single corrugate, top and base protection sheet	50µ LDPE monolayer film tape sealed bag	100µ LDPE monolayer film Heat sealed bag	120µ PE/PA/PE multilayer film vacuum bag heat sealed
Carton dimensions	Carton weight	Units per carton	Units per sterile barrier	
L 550 – W 360 – H 460 mm	6.92 kg	30	5 (PharmaPack 5)	
	Dimensions L x W x H	Weight (inc. pallet)	Cartons per pallet	Units per pallet
Euro pallet	1200x800x1950mm	129kg	16	480
Industry pallet	1200x1000x1920mm	159kg	20	600