

2.7 20L PET (120000-01)

2.7.1 Product Specification Sheet 20L PET

PharmaTainer specifications 20L PET

20L PET PharmaTainer (singlePack)		Catalogue number 120000-01				
Bottle part number	120000-single	USP Class VI	ADCF	DMF No.		
Bottle material	Cellon resin code 001	Yes	Yes	15205		
Nominal volume	20.000ml	Tolerance		Shape (footprint)		
Brim capacity	23.120ml					
Bottle dimensions w/o cap	H500 x W256 x D256				+/- 3.1mm	Square
Bottle weight w/o cap	842.0g				+/- 3%	
Bottle thread O.D.	63.77mm				+/- 0.3mm	
Bottle neck I.D.	54.5mm	+/- 0.3mm	Accuracy			
Graduations	2000ml – 20.000ml	1.000ml intervals	+/- 5%	Endotoxin		
Sterility	C60 γ irradiated X-Ray	SAL 10 ⁻⁶		No detectable Endotoxin Certified by Lot		
Bottle decoration	Graduation scale	Lot number	Expiration date	Data matrix		
Bottle colour	Clear					
Traceability	Lot number	Data matrix	Serial number			
Closure part number	300070C1-99	USP Class VI / EP	ADCF	DMF No.		
Closure material	Cellon resin code 003	Yes	Yes	1646		
		Tolerance				
Closure weight	23.6g	+/-2g				
Closure dimensions	H 28.3mm O.D 70.0mm	+/- 0.3mm				
Closure colour	Natural					
Packaging materials	Carton	Level 3 Carton liner	Level 2	Level 1 Sterile barrier		
	4mm single corrugate, top and base protection sheet	50µ LDPE monolayer film tape sealed bag	100µ LDPE monolayer film Heat sealed bag	120µ PE/PA/PE multilayer film vaccum bag heat sealed		
Carton dimensions	Carton weight	Units per carton	Units per sterile barrier			
L 590 – W 330 – H 570 mm	3.34kg	2	1 (Single Pack)			
	Dimensions L x W x H	Weight (inc. pallet)	Cartons per pallet	Units per pallet		
Euro pallet	1200x800x2190mm	65kg	14	28		
Industry pallet	1200x1000x1840mm	80kg	18	36		