

2.5 5L PET (105000-01)

2.5.1 Product Specification Sheet 5L PET

PharmaTainer specifications 5L PET

5L PET PharmaTainer (singlePack)			Catalogue number 105000-01		
Bottle part number	105000-single	USP Class VI	ADCF	DMF No.	
Bottle material	Cellon resin code 001	Yes	Yes	15205	
Nominal volume	5.000ml	Tolerance	Accuracy	Shape (footprint) Square	
Brim capacity	6.000ml				
Bottle dimensions w/o cap	H322 x W168 x 168 mm				+/- 1.8mm
Bottle weight w/o cap	562.0				+/- 3%
Bottle thread O.D.	63.77mm				+/- 0.3mm
Bottle neck I.D.	54.5mm				+/- 0.3mm
Graduations	500ml – 5.000ml	500ml intervals	+/- 5%	Endotoxin	
Sterility	C60 γ irradiated X-Ray	SAL 10 ⁻⁶		No detectable Endotoxin Certified by Lot	
Bottle decoration	Graduation scale	Lot number	Expiration date	Data matrix	
Bottle colour	Clear				
Traceability	Lot number	Data matrix	Serial number		
Closure part number	300070C1-99		USP Class VI / EP	ADCF	DMF No.
Closure material	Cellon resin code 003		Yes	Yes	1646
			Tolerance		
Closure weight	23.6g		+/-2g		
Closure dimensions	H 28.3mm O.D 70.0mm		+/- 0.3mm		
Closure colour	Natural				
Packaging materials	Carton	Level 3 Carton liner	Level 2	Level 1 Sterile barrier	
	4mm single corrugate, top and base protection sheet	50µ LDPE monolayer film tape sealed bag	100µ LDPE monolayer film Heat sealed bag	120µ PE/PA/PE multilayer film vaccum bag heat sealed	
Carton dimensions	Carton weight	Units per carton	Units per sterile barrier		
L 590 – W 400 – H 350 mm	4.00kg	6	1 (Single Pack)		
	Dimensions L x W x H	Weight (inc. pallet)	Cartons per pallet	Units per pallet	
Euro pallet	1200x800x2140mm	107kg	22	132	
Industry pallet	1200x1000x1910mm	116kg	25	150	