



Liveo™ Pharma Molded Assemblies

Silicone Tubing and Assemblies

Custom moulded silicone assemblies designed for ultra-pure transfer applications in biopharmaceutical and pharmaceutical manufacturing processes.

APPLICATIONS

Liveo™ Pharma Moulded Assemblies are designed for fluid transfer applications requiring integrated premium quality silicone tubing, moulded silicone connections and customer specified components (when applicable).

PROPERTIES AND CUSTOMISATION

Properties of Liveo™ Pharma Moulded Assemblies are dependent upon several factors, including:

- Tubing type: Liveo™ Pharma-50 Tubing, Liveo™ Pharma-65 Tubing, Liveo™ Pharma-80 Tubing, and Liveo™ Pharma Advanced Pump Tubing
- Tubing dimensions: Inside diameter, outside diameter, wall thickness
- Connection type: Sanitary flange end, wye (Y), tee (T), reducer, union, cross, stopper
- Assembly design: Length, number of connections, connection type
- Non-silicone components such as thermoplastic connectors, filters, pinch clamps, and other customer-specified components.

Liveo™ Pharma Moulded Assemblies are available in a range of make-to-order configurations.

Please contact your local DuPont sales representative before writing specifications on this product.



Liveo™ Pharma Molded Assemblies

Silicone Tubing and Assemblies

PACKAGING

Standard packaging includes doublebagging of each assembly in sealed polyethylene bags. Customised packaging options are available; consult your DuPont representative.

Storage and stability

Shelf life

36 months

Characteristics

Food contact

FDA 21 CFR

Additional information

How to use

Information from “Chemical Resistance Guide For Elastomers IV”, published by Compass Publications, La Jolla, California, can be used as a reference. It is however recommended that the DuPont product be evaluated in the real conditions of use.

LIMITATIONS

This product is neither tested nor approved for any hospital or patient care use such as for temporary insertion or any in vivo procedures. This product is not to be used in human implantation, or human contraceptive, reproductive, obstetrical or gynecological applications. The user shall hold DuPont harmless from any and all damages resulting from use of this product. It is the sole responsibility of the user to determine the safety and efficacy of this product for any specific use. This product is not tested for specific pharmaceutical or medical device use(s). Should you wish to use this product in a specific pharmaceutical or medical device application, please contact DuPont to discuss such potential use. It remains the User's responsibility to ensure the safety, efficacy and legal and regulatory compliance in each relevant jurisdiction (including targeted geographic regions of manufacture and supply) of these materials for its intended uses. DuPont makes no representation concerning the suitability of these products for any particular medical or pharmaceutical application. Under no circumstances should these materials be considered for implantation into the human body for periods that exceed 30 days in duration.

Safety warning



Liveo™ Pharma Molded Assemblies

Silicone Tubing and Assemblies

Before handling, please consult the corresponding material safety data sheets. These are available upon request or can be downloaded from our DuPont website : www.DuPont.com

Chemical Media Resistance

Sterilisation methods

- ✓ Autoclave Steam, 30min at 120°C
- ✓ Gamma Radiation, 50 kGy

Symbols used:

- ✓ possibly resistant
Defined as: Supplier has sufficient indication that contact with chemical can be potentially accepted under the intended use conditions and expected service life. Criteria for assessment have to be indicated (e.g. surface aspect, volume change, property change).
- ✗ not recommended - see explanation
Defined as: Not recommended for general use. However, short-term exposure under certain restricted conditions could be acceptable (e.g. fast cleaning with thorough rinsing, spills, wiping, vapor exposure).

DuPont

The information set forth herein is furnished free of charge, is based on technical data that DuPont believes to be reliable, and represents typical values that fall within the normal range of properties. This information relates only to the specific material designated and may not be valid for such material used in combination with other materials or in other processes. It is intended for use by persons having technical skill, at their own discretion and risk. This information should not be used to establish specification limits nor used alone as the basis of design. Handling precaution information is given with the understanding that those using it will satisfy themselves that their particular conditions of use present no health or safety hazards and comply with applicable law. Since conditions of product use and disposal are outside our control, we make no warranties, express or implied, and assume no liability in connection with any use of this information. As with any product, evaluation under end-use conditions prior to specification is essential. Nothing herein is to be taken as a license to operate or a recommendation to infringe on patents.

CAUTION: Do not use DuPont materials in medical applications involving implantation in the human body or contact with internal body fluids or tissues unless the material has been provided from DuPont under a written contract or other acknowledgement that is consistent with the DuPont policy regarding medical applications and expressly acknowledges the contemplated use. For further information, please contact your DuPont representative.

DuPont's sole warranty is that our products will meet our standard sales specifications in effect at the time of shipment. Your exclusive remedy for breach of such warranty is limited to refund of purchase price or replacement of any product shown to be other than as warranted. TO THE FULLEST EXTENT PERMITTED BY APPLICABLE LAW, DUPONT SPECIFICALLY DISCLAIMS ANY OTHER EXPRESS OR IMPLIED WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE, MERCHANTABILITY, OR NON-INFRINGEMENT. DUPONT DISCLAIMS LIABILITY FOR ANY SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES.