



# Clean Image Cleanroom Notebooks

Ideal for use in cleanroom environments.

## Specifications:

Clean Image notebooks are the most advanced line of documentation available for the most critical applications. Clean Image notebooks utilize a Binding Technology that is low in particles and has low particle shed. Ideal for hand writing in a controlled environment. Clean Image notebooks are made of 22lb white college ruled cleanroom paper.



## Anion:

Element	Rating
Chloride (Cl)	0.264 µg/cm <sup>2</sup>
Nitrite (NO <sub>2</sub> )	<0.020 µg/cm <sup>2</sup>
Phosphate (PO <sub>4</sub> )	<0.001 µg/cm <sup>2</sup>
Bromide (Br <sup>-</sup> )	<0.020 µg/cm <sup>2</sup>
Sulfate (SO <sub>4</sub> )	0.272µg/cm <sup>2</sup>

## Cation:

Element	Rating
Lithium (Li)	<0.0015 µg/cm <sup>2</sup>
Sodium (Na)	3.020 µg/cm <sup>2</sup>
Ammonium (NH <sub>4</sub> )	0.03 µg/cm <sup>2</sup>
Potassium (K)	0.15 µg/cm <sup>2</sup>

## College Ruled Paper Specifications

Size	Color	Weight (gsm)	Basis Weight (G/M <sup>2</sup> )	Caliper (micron)	Tear Strength (MD/CD)	Catalog Number	Packaging (Sheets/Notebooks)
8.5"x11	White	22 lbs (80)	78-83	94-104	74/75	PNB CI 8.5x11 2210-L	50sht/10nb, cs
5.5" x 8.5"	White	22 lbs (80)	78-83	94-104	74/75	PNB CI 5.5x8 2220-L	50sht/20nb, cs
3" x 4"	White	22 lbs (80)	78-83	94-104	74/75	PNB CI 3x4 2240-L	50sht/40nb, cs

\*Tensile strength determined by the Elmendorf tear test.

PURUS

THE AMERICAS | +1 760-775-4500 | info@purusint.com

www.purusint.com

PURUS, EnviroTack™, ® or TM are owned by Purus International, Inc. © 2022 Purus International, Inc. All Rights Reserved.



SCAN HERE TO LEARN MORE