

Critical Cover® Coverall® & Combo® Face Shields/Combo Masks & Shield



CRITICAL COVER® COVERALL® FULL & HALF FACE SHIELDS *LIGHT SPLASH PROTECTION*

FEATURES & BENEFITS

When you need total facial splash protection, choose the Coverall® face shield. The face shield's wrap-around design provides the wearer with great peripheral vision while acting as a protective shield against non-hazardous liquid splash and light particles in a Lab or Controlled Environment. The Coverall® face shield is not ANSI Z87 rated for eye protection. You can wear the lightweight shield comfortably for long periods and it fits easily over prescription glasses or protective eyewear. A double-sided anti-fog coating reduces the possibility of fogging in high humidity environments. For safe working conditions, the face shield is low distorting and 100% optically clear while a soft polyurethane foam headband securely positions the device on the wearer.

Available in two styles: full face shield and half face shield. Also available with an adjustable Velcro® band or elastic Comfort headband. Individually packaged.

Coverall Full Face Shields *Styles/Sizes*

Part No.	Description
2802	Full Face Shield w/Velcro Adjustable Band
2807	Full Face Shield w/Comfort Band

100/case (4 boxes of 25)

Coverall Half Face Shields *Styles/Sizes*

Part No.	Description
2803	Half Face Shield w/Velcro Adjustable Band
2808	Half Face Shield w/Comfort Band

100/case (4 boxes of 25)



CRITICAL COVER® COMBO® MASK WITH SHIELD *OPTICALLY CLEAR, REMOVABLE SHIELD*

FEATURES & BENEFITS

The Combo® face mask and face shield is fluid resistant with an optically clear shield. The face shield is attached to the mask with snaps so that it can be easily removed where applicable. The face shield is 7 millimeters thick with an anti-fog coating, low glare and is distortion free. ASTM F2100 Compliant Level 3 mask.

Combo Mask with Shield, Ties *Styles/Sizes*

Part No.	Color
5194-T	Blue

100/case (4 boxes of 25)

Combo Mask with Shield, Ear Loop *Styles/Sizes*

Part No.	Color
5194-E	Blue

100/case (4 boxes of 25)