

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

WHITE POLYPROPYLENE MUFF COVERS



The use of Keystone Cleanroom Product's Polypropylene Muff Covers are an efficient and consistent means of protecting critical surfaces from external sources of particulate contamination after cleaning. The breathable material allows residual moisture from the cleaning process to evaporate.

- Conveniently covers flasks, beakers, containers, media bottles, valves, pipes and hoses
- Prevents dust and other particulate from entering orifices of any size
- Breathable
- Manufactured and packaged in a ISO Certified Cleanroom
- Eliminates Autoclave Tape Residue
- Latex-and talc free rubber elastic
- Lint-Free

PHYSICAL PROPERTIES

Material –100% Non-woven polypropylene Fabric,
CAS# 9003-07-4

Color	White
Basis Weight	1.25 oz. sq. yd

STRENGTH PROPERTIES – (MD machine direction; CD cross direction.)

Elongation (MD)	69.2%
Elongation (CD)	68.4%
Tensile (MD)	99.0 N / 5cm
Tensile (CD)	60.7 N / 5cm

THERMAL PROPERTIES

Melting point is > 225°F.

BARRIER PROPERTIES

Vapor Pressure	negligible mm HG/70F
----------------	----------------------

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

WHITE POLYPROPYLENE MUFF COVERS



KEYSTONE PART NUMBER	FITS OPENING (DIAMETER)
MFNWHD04	1" - 1 3/4"
MFNWHD05	1 3/4" - 2 1/2"
MFNWHD06	2" - 3 3/4"
MFNWHD07	2 1/2" to 4"
MFNWHD08	3 1/2" to 4 3/4"
MFNWHD09	4" to 5 1/2"
MFNWHD12	5 1/2" to 8"
MFNWHD16	8" to 12"
MFNWHD18	10" to 14"